INSTRUCTIONS FOR USE

Ambu® aBox™ Duodeno










For use by trained clinicians/physicians only.

For in-hospital use.

For use with Ambu® aScope™ Duodeno.

Ambu

CE 2797

Ambu is a registered trademark and aBox™ Duodeno is a trademark of Ambu A/S.

Contents

Page

English (Instruction for use).....	4-27
Български (Инструкция за употреба).....	28-55
Čeština (Návod k použití).....	56-81
Dansk (Brugervejledning).....	82-105
Deutsch (Bedienungsanleitung).....	106-132
Ελληνικά (Οδηγίες χρήσης).....	133-160
Español (Instrucciones de uso).....	161-187
Eesti (Kasutusjuhend).....	188-212
Suomi (Käyttöohje).....	213-236
Français (Mode d'emploi).....	237-261
Hrvatski (Upute za upotrebu).....	262-287
Magyar (Használati útmutató).....	288-312
Italiano (Istruzioni per l'uso).....	313-337
日本語 (使用説明書).....	338-360
Lietuviškai (Naudojimo instrukcija).....	361-383
Latviski (Lietošanas instrukcija).....	384-408
Nederlands (Gebruiksaanwijzing).....	409-433
Norsk (Bruksanvisning).....	434-456
Polski (Instrukcje użytkowania).....	457-482
Português (Instruções de utilização).....	483-508
Română (Instrucțiuni de utilizare).....	509-533
Русский (Инструкция по применению).....	534-562
Slovenčina (Návod na použitie).....	563-587
Slovenščina (Navodila za uporabo).....	588-612
Svenska (Bruksanvisning).....	613-635
Türkçe (Kullanma talimatları).....	636-659
中文 (使用说明).....	660-679

1. Important Information – Please Read Before Use!

Read this *Instructions for Use* (IFU) before operating and keep for future reference. Failure to read and thoroughly understand the information presented in this IFU, as well as those developed for ancillary endoscopic equipment and accessories, may result in serious injury to the patient and/or user. Furthermore, failure to follow the instructions in this IFU may result in damage to, and/or malfunction of, the equipment.

This IFU describes the recommended procedures for inspecting and preparing the equipment prior to its use. It does not describe how an actual procedure is to be performed, nor does it attempt to teach the beginner the proper technique or any medical aspects regarding the use of the equipment. It is the responsibility of each medical facility to ensure that only appropriately trained personnel, who are competent and knowledgeable about the endoscopic equipment, antimicrobial agents/processes and hospital infection control protocol are involved in the use, handling and the care of these medical devices. The IFU may be updated without further notice. Copies of the current version are available upon request.

Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (e.g. IEC 60950 or IEC 62368 for data processing equipment). Furthermore, all configurations shall comply with the requirements for medical grade electrical systems (see clause 16 of the latest valid version of IEC 60601-1). Anybody connecting additional equipment to medical grade electrical equipment which configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. If in doubt, consult your local representative or the technical service department.

US federal law restricts this device for sale only by, or on the order of, a physician.

1.1. Intended Use/Indication for Use

The aBox™ Duodeno is designed to be used with the aScope™ Duodeno, and other ancillary equipment (e.g. medical grade video monitor) for endoscopy and endoscopic surgery within the duodenum.

NOTE: Do not use this device for any purpose other than its intended use. Select the endoscope to be used according to the objective of the intended procedure based on the full understanding of the endoscope's specifications and functionality as described in this IFU.

1.2. Target User Groups and User Qualifications

The device is intended to be used by physicians trained on ERCP procedures. If there are official standards for user qualifications to perform endoscopy and endoscopic treatment that are defined by the hospital's medical administrators or other official institutions, such as academic societies on endoscopy, follow those standards. If there are no official qualification standards, the operator of this device must be a physician approved by the medical safety manager of the hospital or person in charge of the department (e.g. department of internal medicine, etc.).

The physician should be capable of safely performing the planned endoscopy and endoscopic treatment following guidelines set by the academic societies on endoscopy, etc., and considering the difficulty of endoscopy and endoscopic treatment. This manual does not explain or discuss endoscopic procedures.

1.3. Contraindications

Contraindications depend on the endoscope used and the endoscopic procedure. For detailed information regarding contraindications please refer to the IFU of the aScope™ Duodeno.

1.4. Installation and Maintenance

The medical devices described in this IFU must be tested/inspected in accordance with national regulations during installation and regular inspection. The medical device does not require regular maintenance.

1.5. Warnings and Cautions

Failure to observe these warnings and cautions may result in patient injury or damage to the equipment. The manufacturer is not responsible for any damage to the system or patient injury resulting from incorrect use.

WARNINGS

Indicate a potentially hazardous situation which, if not avoided, could result in death or serious injury.

Preparation for Use

- Do not connect non-medical grade devices to the aBox™ Duodeno. Connection of devices that are not medical grade (IEC 60601 conform) could have negative impact on the safety of the system. Connect only medical grade equipment to the aBox™ Duodeno.
- Never use the aBox™ Duodeno if an abnormality is suspected. Damage or irregularity in the device may compromise patient and/or user safety and may result in more severe equipment damage.
- Keep fluids away from all electrical equipment. If fluids are spilled on or into the unit, stop operation of the aBox™ Duodeno immediately and contact Ambu. Do not prepare, inspect or use the aBox™ Duodeno with wet hands.
- In case of device failure or malfunction, always keep another aBox™ Duodeno in the room ready for use.
- Never insert or spray anything into the ventilation grills of the aBox™ Duodeno. It can cause an electric shock and/or fire.
- Always set the minimum required brightness. The brightness of the image on a medical grade video monitor may differ from the actual brightness at the distal end of an endoscope. Although the illumination light emitted from the endoscope's distal end is required for endoscopic observation and treatment, it may also cause alteration of living tissues such as protein denaturation of liver tissue and perforation of the intestines by inappropriate use.
- Do not leave the endoscope illuminated before and after examination. Otherwise, the LEDs could be redirected to the opened eyes of the sedated patient and can cause retina burns. This product may interfere with other medical grade electronic equipment used in combination with it.
- Before use, refer to the Appendix 1 to confirm the compatibility of this device with all equipment to be used.
- Do not use this product in any place where it may be subject to strong electromagnetic radiation (for example in the vicinity of a microwave therapeutic device, MRI), wireless set, short-wave therapeutic device, cellular/portable phone, etc.). This may impair the performance of the product.
- If the endoscopic image dims during use, blood, mucus or debris may adhere to the light guide on the distal end of the endoscope. Try to clear the LEDs by rinsing. If image is still dimmed, carefully withdraw the endoscope from the patient and remove the blood or mucus in order to obtain optimum illumination and to ensure the safety of the examination. If you continue to use the endoscope in such a condition, the distal end temperature may rise and cause mucosal burns. It may also cause patient and/or operator injury.
- To display observation images, connect the output terminal of the aBox™ Duodeno directly to the monitor. Do not make the connection via any supplementary image- or video communication equipment device not provided by Ambu. Images may disappear during observation depending on the condition of such supplementary equipment.

Inspection of the aBox™ Duodeno

- Do not connect the power plug to the 2-pole power circuit with a 3-pole to 2-pole adapter. To avoid the risk of electric shock, this equipment must only be connected to a mains supply with protective earth.

- Do not use the aBox™ Duodeno if not inspected as instructed. Inspect other equipment to be used with the aBox™ Duodeno as instructed in their respective instruction manuals. Should any irregularity be observed, do not use the aBox™ Duodeno and see section 12. Troubleshooting. If the irregularity is still observed after consulting section 12, contact Ambu. Damage or irregularity may compromise patient or user safety and may result in more severe equipment damage.
- Do not use this device when the live image cannot be observed. Otherwise, patient injury may occur.

Connection to the AC Mains Power Supply

- Always keep the power plug dry. A wet power plug may cause electric shocks.
- Confirm that the hospital-grade wall mains outlet to which this device is connected has adequate electrical capacity that is larger than the total power consumption of all connected equipment. If the capacity is insufficient, fire can result or circuit breaker may trip and turn OFF this device and all other equipment connected to the same power circuit.
- Be sure to connect the power plug securely to prevent erroneous unplugging during use. Otherwise, the equipment will not function.
- If combinations of equipment other than those shown below are used, the full responsibility should be assumed by the medical treatment facility. Such combinations do not only allow the equipment to manifest their full functionality but may also imperil the safety of the patient and medical personnel. In addition, the endurance of the video system center and ancillary equipment is not guaranteed. Troubles caused in this case are not covered by free-of-charge repair. Be sure to use the equipment in one of the recommended combinations.

Operation of the aBox™ Duodeno

- To guard against dangerous chemicals and potentially infectious material during the procedure and danger of unintentional diathermy burns, wear personal protective equipment, such as eyewear, face mask, moisture-resistant clothing and chemical and electrical-resistant gloves that fit properly and are long enough so that your skin is not exposed. Please note that a new pair of gloves is required prior to each procedure.
- Never use the aBox™ Duodeno if an abnormality is suspected. Damage or irregularity in the device may compromise patient or user safety and may result in more severe equipment damage.
- If any other abnormality occurs or is suspected, immediately stop using the equipment, turn OFF all equipment, and gently withdraw the endoscope from the patient as described in the endoscope's instruction manual. Then, refer to the instructions in section 12. Troubleshooting. If the problems cannot be resolved by the remedial action described in section 12, do not use the equipment and contact Ambu.
- Clean the device and change gloves prior to touching the device and between cases. Not changing gloves could lead to cross-contamination.

Ancillary Equipment

- When using spray-type medical agents such as lubricant, anesthetic, or alcohol, use them away from the aBox™ Duodeno so that the medical agents do not come into contact with the aBox™ Duodeno. Medical agents might leak into the video system through the ventilation grills and may cause equipment damage.
- Do not use non-compatible electrosurgical equipment with this device. Interference on the monitor or loss of the endoscopic image can occur.
- Do not use a humidifier near the video system center as dew condensation might occur and may cause equipment failure.
- When recording images, be sure to record the images together with patient data. Otherwise, distinction between different observations may become difficult.

Fuse Replacement

- Never use a fuse other than the fuse model designated by Ambu. Otherwise, malfunction or failure of the aBox Duodeno may cause a fire or electric shock hazard.
- Be sure to turn the aBox™ Duodeno OFF and unplug the power cable before removing the fuse from the aBox™ Duodeno. Otherwise, fire or electric shock may result.
- If the power fails to come on after replacing the fuses, unplug the power cable immediately from the AC mains power inlet and then contact Ambu. Otherwise, electric shock may result.
- Do not position the device above the patient. If the front is destroyed due to strong external impacts, splintering may occur. Falling splinters could cause injury to the patient.
- Insert the fuse box into this device until it clicks into position. If the fuse box is inserted incompletely, the power may fail to come ON or a power failure may occur during operation.

Care and Storage

- After wiping with a piece of moistened gauze, dry the aBox™ Duodeno thoroughly before using it again. If it is used while still wet, there is the risk of an electric shock.
- Do not use the device without thoroughly understanding the IFU. If the system is not properly prepared before each use, equipment damage, patient and operator injury and/or fire can occur.

Potential of Fire

- The system is not protected against fire and explosion. When using the device in areas with flammable or explosive gases or in areas with oxygen-enriched air, fire or explosions may occur in the unit. Do not operate the device in areas with flammable or explosive gases or gas mixtures. Do not operate the device in an oxygen-enriched environment.

Installation, Repair and Maintenance

- The device does not contain any parts that can be repaired by the user. Any disassembling, change or attempt to repair can lead to injury of patient or user and damage of the system. Installation, repairs and maintenance must be performed by Ambu staff or Ambu approved staff. Please read section 12 for more information regarding troubleshooting.

CAUTIONS

Indicate a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices or potential equipment damage.

General

- Do not use a pointed or hard object to press the buttons on the front panel. This may damage the buttons.
- Do not apply excessive force to the aBox™ Duodeno and/or other devices connected. Otherwise, damage and/or malfunction can occur.
- Clean and vacuum dust the ventilation grills of the aBox™ Duodeno using a vacuum cleaner, when necessary. Keep the ventilation grill clear. Otherwise, the aBox™ Duodeno may break down and get damaged from overheating.
- Be sure that this device is not used adjacent to or stacked with other equipment (other than the components of this system) to avoid electromagnetic interference.
- Electromagnetic interference may occur to this device when it is placed near equipment marked with the following symbol (Ⓜ) or other portable and mobile RF communications equipment such as cellular phones. If radio interference occurs, mitigation measures may be necessary, such as reorienting or relocating this device or shielding the location.

- Do not place any object on the top of the aBox™ Duodeno. Otherwise, equipment deformation and damage can result.
- Place the aBox™ Duodeno on a stable, level surface. Otherwise, the aBox™ Duodeno may topple down or drop, and user or patient injury may occur or equipment damage can result.
- If a trolley other than the mobile workstation is used, confirm that the trolley can withstand the weight of the equipment installed on it.
- US federal law restricts this device for sale only by, or on the order of, a physician.

Care and Storage

- Do not clean the power cable socket, the connections and the AC mains power inlet. Cleaning can deform or corrode the contacts, which could damage the aBox™ Duodeno.
- Do not store the device in a location exposed to direct sunlight, X-rays, radio activity or strong electromagnetic radiation (e.g. near microwave medical treatment equipment, short-wave medical treatment equipment, MRI equipment, radio or mobile phones). Otherwise, damage to the aBox™ Duodeno can occur.
- When disposing of this device or any of its components (such as fuses), follow all applicable national and local laws and guidelines.
- Turn OFF all ancillary equipment before connecting to the aBox™ Duodeno and use appropriate cables only. Otherwise, equipment damage or malfunction can result.
- Close the connector cover before you clean the aBox™ Duodeno. Do not open the cover while cleaning the aBox™ Duodeno, otherwise fluid can penetrate the connector and damage the device.
- Do not autoclave or sterilize the aBox™ Duodeno. This could lead to a damage of the device.
- The cables should not be sharply bent, pulled, twisted or crushed. Cable damage can result.
- Never apply excessive force to connectors. This could damage the connectors.
- Use this device only under the conditions described in "Transportation, storage and operation environment Specifications" in section 10.2. Otherwise, improper performance, compromised safety and/or equipment damage may result.

Accessories

- Do not attach anything other than the bottle holder to the aBox™ Duodeno. Otherwise, equipment damage or deformation may result.
- Only place a bottle of sterile water into the bottle holder. Otherwise, the bottle holder can be damaged.
- When the device is used with energized endoscopic devices, leakage current may be additive. Use only endoscopic devices of type BF or CF. Check the compatibility of the accessory/endoscopic device before use regarding any criteria for safe use.

1.6. Adverse Events




No known adverse events for the aBox Duodeno, please refer to the aScope Duodeno IFU for ERCP related adverse events.

Potential adverse events in relation to Ambu Duodeno System (not exhaustive): Infection/ inflammation (including post-ERCP pancreatitis (PEP), cholangitis, cholecystitis, endocarditis and sepsis), bleeding, perforation, thermal injuries, stent related adverse events, cardiopulmonary adverse events, air embolism, anaesthesia related adverse events, nausea, sore throat, abdominal pain and discomfort.

2. System Description


2.1. System Parts

The aBox™ Duodeno is reusable. No modification of this equipment is allowed. The aBox™ Duodeno is shipped with one power cable that supplies the power necessary to run the aBox™ Duodeno, a bottle holder for the sterile water bottle and a potential equalization cable (POAG).

Ambu® aBox™ Duodeno Reusable Device	Part Number
	485001000US (for US market) 485001000 (for non US markets)
Bottle Holder	
	
Power Cable	
	
POAG Cable	
	

The aBox™ Duodeno is not available in all countries. Please contact your local sales office.


2.2. Compatible Device (Application Part)

Ambu® aScope™ Duodeno Sterile and single-use device	Part Number	Distal End Outer Diameter \varnothing	Working channel Inner Diameter \varnothing
	482001000	41.2 Fr (13.7 mm)	12.6 Fr (4.2 mm)

The aScope™ Duodeno is not available in all countries. For detailed information please contact your local sales office.

2.3. aBox™ Duodeno Description and Function

The aBox™ Duodeno is the console that is necessary for processing the endoscope camera's video image, the remote switch signals and output video and recorder data. The aBox™ Duodeno is designed to be used with the aScope™ Duodeno.

No.	Part	Function
1	Bottle holder (sterile water)	Retainer for sterile water
2	Rinsing pump (peristaltic pump)	Lens rinsing
3	Connector snap	Bracket for endoscope connector
4	Endoscope connector cover	Protection cover
5	Endoscope connector socket	Electrical connection between aScope™ Duodeno and aBox™ Duodeno
6	Illumination button	Press the button for illumination
7	Control panel	Buttons for navigation through the settings and information menu of the aBox™ Duodeno.
8	Power button with power indicator	Press the button for power ON before procedure and power OFF after procedure. The button will light up green when turned ON
9	Main fuse	Device protection
10	Main power switch	Switch button for power ON and power OFF
11	AC power inlet	AC power cable socket
12	POAG socket	Socket for potential equalization
13	System ventilation	Ventilation
14	HD-SDI connections	Video output
15	DVI connection	Video output
16	Stereo jack 3.5 mm	Trigger output for video and image capture
17	D-SUB 9P	Trigger output for video and image capture
18	Service panel	Access for service technicians

3. Explanation of Symbols Used

3.1. Symbols

Symbols	Indication	Symbols	Indication
	Consult <i>Instructions for Use</i>		NRTL Symbol
	Protect the packaged product from moisture		Medical Device
	Waste Bin symbol, indicating that waste must be collected according to local regulation and collection schemes for disposal of electronic and electrical waste (WEEE)		Protection against electric shock - Type BF, safety class IEC60601-1
			Refer to <i>Instructions for Use</i>
IP 21	Ingress protection	Rx Only	Prescription Device
	Serial number (consisting of number and year of manufacturing)		Reference Number
	Ground Connection		Potential Equalization
	Legal Manufacturer		Manufacturing Site
	CE Mark		Do not touch moving parts

4. Ancillary Equipment and Accessories

The following ancillary equipment must be connected for the system to become functional and for recording patient data. Inspect the following equipment as described in their respective instruction manuals.



- Bottle holder (part of the system).
- Sterile water for rinsing equal to or greater than 1000 ml in volume.

It is recommended that the bottle of sterile water is placed in the designed bottle holder on the left side of the aBox™ Duodeno. For detailed information regarding attaching the bottle holder to the aBox™ Duodeno please refer to the section 4.1 below.

- Medical Grade monitor with resolution of at least 1920x1080 and a monitor size of at least 27" with DVI or HD-SDI input(s). Recommended color space is sRGB.
- Image Capture Report and/or writing workstation.
- Vacuum source of at least -7 psi (-50 kPa) with suction system.
- Insufflation source approved for use in GI endoscopic procedures.

4.1. Attaching the Bottle Holder to the aBox™ Duodeno

The aBox™ Duodeno is shipped with one bottle holder that must be attached to the left side of the device. Follow the illustrations below to properly attach the bottle holder.

5. Preparation and Inspection for Use

Numbers in gray circles below refer to illustrations on page 2. Inspect all ancillary equipment to be used with this device as instructed in their respective *Instructions for Use*. Should any irregularity be observed after inspection, follow the instructions as described in section 12. Troubleshooting. If this device malfunctions, do not use it. Contact your Ambu sales representative for further assistance.

- Inspect the content of the aBox™ Duodeno. Match all items in the package with the components shown in the device description in section 2. **1**
- If the device is damaged, a component is missing or you have any questions, do not use the device, immediately contact Ambu.
- It is recommended that the aBox™ Duodeno is placed on level surfaces with the capability of movement (i.e. mobile cart, medical equipment booms) such that the system can be moved to the most advantageous position for any given patient and/or as required to satisfactorily perform the intended patient procedure. All such carts or booms should be designed for this purpose and rated for the weight requirements necessary (refer to the technical data section 10, as well as possessing a locking mechanism to prevent inadvertent rolling or movement during a procedure. **2**)
- Attach the Bottle Holder to the aBox™ Duodeno as described in section 4.1. **2**

- One (1) hospital grade power supply cable is provided with the aBox™ Duodeno and is necessary to supply power from the electrical mains to the aBox™ Duodeno (refer to the technical data section 10 for electrical ratings and other applicable information). The power supply cable is not an integral part of the aBox™ Duodeno. Connect the power supply cable to the AC Mains Connection and to a grounded power source socket. **3**
- This medical device may be connected to a network of medical grade devices. Use the potential equalization conductor as determined necessary by your facilities Biomedical/Clinical/Technical Engineering staff. The potential equalization conductor (is easily identifiable as the green cable with a yellow line running along its length) serves as a conductor for possible differences in ground potentials between network components which can result in leakage current that can flow to the patient and is potentially dangerous. The potential equalization conductor works to remove this hazard. **3**
- Connect the aBox™ Duodeno to at least one medical grade monitor. It is recommended to use HD-SDI 1 for the main monitor and to choose a medical grade monitor with a full HD resolution. Please refer to Appendix 1 for specific connections details. **3**
- An additional medical grade monitor or a medical grade recorder can be connected to HD-SDI 2 or DVI. **3**
- Use the "TRIG 1, TRIG 2 and / or TRIG 3" output for the remote triggering signal cable to the medical grade recorder whenever a recorder is connected. For detailed information regarding connections please refer to Appendix 1. **3**

5.1. Powering up and Starting the aBox™ Duodeno

Once all components as described above have been connected, the aBox Duodeno can be powered up.

- Switch the aBox™ Duodeno ON by switching on the main switch on the back of the unit and then by pressing the power button on the front (right) side of the unit. The power button will light up green when the aBox™ Duodeno is powered up. **4a**
- The aBox™ Duodeno indicates when it is ready by displaying the information message on the main monitor screen: "please connect endoscope". Once the endoscope is connected, the system will confirm by displaying the message "endoscope connected" on the main monitor screen.

5.2. Preparing and Connecting the Ambu® aScope™ Duodeno

Please refer to the IFU of the aScope™ Duodeno. **5**


5.3. The Main Monitor Screen

After completing the steps described in sections 4. and 5., the system will perform an internal system check and display the status on the main monitor screen.

You will now have to confirm that a new bottle of sterile water has been prepared prior to the procedure and you will be asked to connect the CO₂ and the vacuum source as indicated. Once the connections have been made, you can perform the functional check of the endoscope. To confirm, please press the  button on the control panel of the aBox™ Duodeno. After confirmation, the aBox™ Duodeno will display the live image on the main monitor screen. Press the illumination  button for light. **6**

5.4. Control Panel Navigation

The settings and information menu can be displayed by pressing any button on the control panel of the aBox™ Duodeno.

Explanation of control panel buttons

Button	Name	Function
	Illumination button	Turns the aScope™ Duodeno LEDs ON and OFF
	Left button	Navigate to the left. The left button is used to hide informations and settings menu
	Right button	Navigate to the right
	Up button	Navigate upwards
	Down button	Navigate downwards
	OK button	Confirms messages and/or settings



Image settings ①. Open the menu by pressing any button on the control panel, set up and then go back with the **left button** . Press once to exit the setting and twice to close the menu.

Remote switches ②. Open the menu by pressing any button on the control panel and select Remote switches. The available function for each remote switch will be displayed when entering the settings in a pull side menu submenu. Set the desired function for each remote switch and save.

- To exit the information and settings menu, use the left button (◀). Press once to exit the setting and twice to close the menu.
- The remote switches standard configuration is: 1 for image capture, 2 and 3 have no preset function.

Language (3). Notice the language setting is only active in the idle mode. If a scope is connected you can't change the language.

Open the menu by pressing any button on the control panel, navigate with **down button** (▼) to **Language** (3) and enter in with the **right button** (▶). Choose your language by using the left or right buttons. The language will be changed directly. Confirm with (OK) button.

Shutdown (4). After the final procedure of the day, select the shutdown button on the settings and information menu. After selecting to shutdown, press the shutdown button for ending the daily session, press the (OK) button on the control panel to confirm. The system will now begin the process of powering down. Shutdown is completed when the power button on the aBox™ Duodeno is no longer illuminated.

Info button (5). When setting up the aBox™ Duodeno, you will be asked to confirm that a new bottle of sterile water has been prepared pre procedure. The info option will save every message which might pop-up during the procedure, e.g. "Lid of the peristaltic pump is open". Press the (OK) button to confirm and the **left button** (◀) to close the menu.

Standard button (6). Press the Standard button in any Sub Menu to return to standard settings of each Sub Menu.

The system will display messages on the main screen monitor if an irregularity occurs. For detailed information, see section 12. Troubleshooting.

Explanation of the Symbols		
Symbol	Name	Description
	Peristaltic pump	Icon will be displayed when the lid of the peristaltic pump is open
	Warning	Icon appears if a message is displayed and not confirmed
	LEDs ON	The number in this symbol will indicate the level of increased or decreased brightness
	LEDs OFF	The number in this symbol will indicate the level of increased or decreased brightness

Explanation of the Symbols		
Symbol	Name	Description
	Image capture	This icon will be displayed when an image has been captured
	Video recording	This icon will be displayed when video recording begins
	Zoom	This icon will be displayed when the zoom function is ON
	Rinsing	This icon will be displayed when the rinsing function is activated

6. Ending a Procedure and Shutting Down the System

6.1. Concluding a Patient Procedure

To conclude a patient procedure, remove the endoscope connector plug from the connector socket of the aBox™ Duodeno. For detailed information and disposal of the endoscope, please refer to the IFU of the aScope™ Duodeno. **7**

6.2. System Shutdown

After the final procedure of the day, push the power button for 3 seconds. The system will now begin the process of powering down. Shutdown is completed when the power button on the aBox™ Duodeno is no illuminated.

Notice that incorrect shut-down of the aBox™ Duodeno may impair its functionality permanently and may require service. **8**


7. Fuse Replacement

Always use the fuses recommended by Ambu. To order new fuses, contact your sales representative. Make sure that only appropriately trained personnel are in charge of the fuse replacement.

Turn the aBox™ Duodeno OFF and disconnect the power cable from the mains outlet.

- Release the fuse by squeezing the lower tab of the fuse box using a screwdriver. **1**
- Pull the fuse box straight out. **2**
- Replace both fuses. **3**
- Insert the fuse box into the aBox™ Duodeno until it clicks into position. **4**


Connect the power cable to a grounded power source socket and turn the aBox™ Duodeno ON and confirm the power output.

8. Peristaltic Pump Replacement

Ensure the aBox™ Duodeno is powered off and that the peristaltic pump housing is completely closed.

- Release the entire peristaltic pump head by pressing the release tab on the lower-right side of the pump housing. ①
- Rotate the peristaltic pump to your left and pull it towards you in one continuous motion. ②
- Align the key of the replacement peristaltic pump assembly with the corresponding notch on the aBox™ Duodeno, then push and rotate the pump assembly to the right until it locks into position. ③
- When properly secured, peristaltic pump head should not rotate in any direction. ④

9. Cleaning of the aBox™ Duodeno External Surfaces

The aBox™ Duodeno should be cleaned before and after each procedure. Clean the aBox™ Duodeno according to good medical practice and using the procedures below: 9

Use a germicidal wipe to remove heavy soil. All blood and other body fluids must be thoroughly cleaned from all surfaces and objects before using germicidal wipe.

(SUPER- SANI-CLOTH® from PDI or equivalent)

The treated surfaces must remain visibly wet for the full length of two (2) minutes. Use additional wipes, if needed, to assure continuous 2 minutes wet contact time. Let the aBox™ Duodeno air dry.

10. Technical Data

All the following reported measurements (e.g. weight, dimensions) are average values. Therefore, small variations may occur, which however have no effect on the performance and safety of the system.

10.1. aBox™ Duodeno Specifications

Power Supply	Voltage	AC 120 V / 230 V
	Frequency	50 / 60 Hz
	Consumption electric power	91 VA
	Fuse rating	2x 5 A H 250 V T
	Fuse size	5 mm x 20 mm
Size of the aBox™ Duodeno	Dimensions	494 (D) x 487 (W) x 145 (H) mm
	Weight	13 kg

Classification (medical electrical equipment)	Type of protection against electric shock	Protection Class I
	IP Classification	IP21

10.2. Transportation, Storage and Operation Environment Specifications

Transportation temperature	-5 °C – 40 °C (23 °F – 104 °F)
Storage temperature	10 °C – 25 °C (50 °F – 77 °F)
Operation temperature	10 °C – 40 °C (50 °F – 104 °F)
Relative humidity transportation and operation	30 – 85 % relative
Relative humidity storage	10 – 43 % relative
Atmospheric pressure	80 – 109 kPa (100 kPa = 1 Bar) 11.6 – 15.8 psi 600 – 818 mm Hg

10.3. Accessories

General information

Connected equipment, especially electrical equipment, must conform to relevant medical standards (medical grade) as described in section 4.

10.4. List of Accessories-Examination

Accessory	Information	Part of the System
Rinsing water	Sterile water available in clinical environment, with at least 1000ml	No
Bottle holder	Holder for bottles of sterile water	Yes
Suction canister	Medical grade vacuum suction canister. Any marketed canister can be used	No

11. Environmental Protection

With regards to European Union Directive 2002/96/EC on waste electrical and electronic equipment (WEEE), all medical waste electrical and electronic equipment (WEEE) should be disposed of and collected separately. This product is Electrical and Electronic Equipment and should be disposed of in accordance with national and local legislation and requirements.

12. Troubleshooting

The following table shows the possible causes and countermeasures against troubles that may occur due to equipment setting errors or deterioration of the **aBox™ Duodeno**. Troubles or failures other than those listed in the following table need repair. As repair performed by persons who are not qualified by Ambu could cause patient or user injury and/or equipment damage, be sure to contact Ambu for repair.

Irregularity Description	Possible Cause	Solution
Power fails to come on	The power switch and/or the power button of the aBox™ Duodeno is set to OFF	Set the power switch to ON
	The power cable is not connected	Connect the power cable with an electrical source as described in section 5
No video image	The aBox™ Duodeno or ancillary equipment is not switched ON	Switch aBox™ Duodeno and ancillary equipment on
	Medical grade monitor not properly connected or defective	Connect the medical grade monitor properly
	aScope™ Duodeno not properly connected or defective	Connect the aScope™ Duodeno correctly or connect a new aScope™ Duodeno
	LEDs not switched ON	Switch the LEDs ON
aBox™ Duodeno does not boot up	aBox™ Duodeno not powered ON	Power the aBox™ Duodeno ON
	Medical grade monitor not properly connected or defective	Shutdown aBox™ Duodeno and connect the medical grade monitor properly. Try to power ON the aBox™ Duodeno after 10 sec again
Power fails to come on	Power cable not connected	Connect the power cable
	Fuse blown	Replace fuse. Please refer to section 7
aBox™ Duodeno will not shut down	aScope™ Duodeno still connected	Please unplug the aScope™ Duodeno from the aBox™ Duodeno and press the power button for 3 secs. to shut down the system
	Power button was pressed for too short a duration	Press the power button for at least 3 sec
	aBox™ Duodeno defective	Switch the aBox™ Duodeno off by using the mains switch on the back panel of the aBox™ Duodeno and contact Ambu for support
Rinsing not possible	Peristaltic pump defective	Replace the peristaltic pump. See section 8

Irregularity Description	Possible Cause	Solution
Control panel does not work	aBox™ Duodeno not powered ON	Power the aBox™ Duodeno ON
	Sterile water bottle not confirmed	Confirm the message on the settings and information menu that a new bottle of sterile water was connected
	aScope™ Duodeno not connected	Connect the aScope™ Duodeno
Endoscopic image is too dark	LEDs are not switched ON	Switch the LEDs on as described in section 5
	Image settings of the aBox™ Duodeno are incorrect (brightness and contrast)	Set the image settings correct as described in section 5.4
	Medical grade monitor settings are improper (brightness and contrast)	Set a proper brightness as described in the medical grade monitor's instruction manual.
	LEDs are running in reduced power mode	Check the sterile water bottle and connect a new bottle of sterile water if necessary or check if the process water tubing is sufficiently inserted in the sterile water bottle
Check the suction system for proper functionality. See section 12. Troubleshooting in the aScope™ Duodeno IFU		
Endoscopic image is too bright	Image settings of the aBox™ Duodeno are incorrect (brightness and contrast)	Set the image settings correct as described in section 5.4
	Monitor settings are improper (brightness and contrast)	Set a proper brightness or contrast as described in the monitor's instruction manual
Color tone of the endoscopic image is unusual	Medical grade monitor settings are improper	Adjust the color setting at the monitor according its IFU, starting from standard color setting (D65)
	Monitor cable is connected incorrectly	Connect the monitor cable properly as described in Section 5
	Medical grade monitor cable is defect	Connect a new power cable to the medical grade monitor
Endoscopic image remains frozen	Freeze button is still set	Press the freeze button to restore the real-time image
Image cannot be stored	No video recorder connected	Connect a video recorder

Standard Messages During Set-up and Procedure

Message	Possible Cause	Solution
Please connect a new bottle of sterile water	Standard request before each examination	Confirm, that a new bottle of sterile water was connected
Lid of the peristaltic pump is open	Lid of the peristaltic pump is open	Close the lid of the peristaltic pump
Please check process water and the connection of the suction pump	Increased temperature of the distal end	Check process water and the connection of the suction pump
Image initializing. Please wait	Standard message after an image drop-out was detected	Wait until the aBox™ Duodeno displays an image on the monitor. If nothing happens, contact Ambu for support

Error Messages

Message	Possible Cause	Solution
Video signal issue identified (E1)	aScope™ Duodeno defect	Connect a new aScope™ Duodeno
	Electrosurgical device set on increased intensity	Reduce the intensity of the electrosurgical device
aBox™ Duodeno defective (E2)	Framer grabber not connected/defective/detected	End the procedure and contact Ambu for support
aBox™ Duodeno defective (E3)	No communication between the controller board and the aBox™ Duodeno	End the procedure and contact Ambu for support
aBox™ Duodeno defective (E4)	Wrong settings in the controller board.	End the procedure and contact Ambu for support
Remote switches and rinsing function access limited.	Limited functionality of the endoscope due to a defect endoscope	If necessary, take a new aScope™ Duodeno for concluding the procedure
aBox™ Duodeno defective (E5). Contact Ambu for support	Hardware and software error	End the procedure and contact Ambu for support

12.1. Returning the aBox™ Duodeno to Ambu

Should it be necessary to return an aBox™ Duodeno to Ambu for investigation, please contact your Ambu representative beforehand for instructions and/or guidance. To prevent infection, it is strictly forbidden to ship contaminated medical devices outside of our guidance. Hence, the aBox Duodeno must be appropriately cleaned and packed on site before shipment to Ambu. Ambu reserves the right to return contaminated medical devices to the sender.

In an event of a serious incident, please inform Ambu and the competent authority.

12.2. Standards Applied

The system conforms to the standard of IEC 60601-1 / ANSI/AAMI ES60601-1, IEC 60601-1-2 and IEC 60601-2-18.

12.3. Electromagnetic Compatibility

General Information Medical electrical equipment is subject to special precautions with respect to EMC and must be installed according to the instructions in the accompanying documentation.

The manufacturer can only guarantee compliance of the equipment if the accessory parts listed in the accompanying documentation are used.

The device is intended for use exclusively by trained medical personnel. This device can cause radio interference or interference with the operation of other equipment in its close vicinity. It may be necessary to take suitable corrective measures, such as readjustment, a different system layout, or shielding.

Special instructions Regulations applying to medical devices require to provide you with the following information. (See the following pages of all tables.

- Table 1 Recommended safety distances.
- Table 2 Electromagnetic compatibility 1.
- Table 3 Electromagnetic compatibility 2.
- Table 4 Electromagnetic transmission.

Table 1 Recommended Safety Distance

Table 1 Recommended Safety Distance			
Recommended safety distance between portable and mobile HF telecommunications systems and the aScope™ Duodeno. The system is intended for operation in an electromagnetic environment in which HF interference is controlled. The user of the system can prevent electromagnetic interference by maintaining a minimum distance between portable and mobile HF telecommunications systems and the System, depending on the output power of the communication unit, as specified below.			
Nominal power, P, of transmitter, measured in watts [W]	Recommended safety distance, d, expressed in meters, based on the nominal transmitter power and transmission frequency		
	150 kHz – 80 MHz $d=3.5/3 \sqrt{P}$	80 MHz – 800 MHz $d=3.5/3 \sqrt{P}$	800 MHz – 2.5 GHz $d=3.5/3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	12	12	23
For transmitter where the maximum nominal power is not specified, the recommended safety distance can be determined by use of the formulas above.			

Table 2 Electromagnetic Compatibility 1

Guidelines and manufacturer's declaration – resistance to electromagnetic interference. The aScope™ Duodeno is intended for operation in the electromagnetic environment specified below. The user of the aScope™ Duodeno should ensure that it is used in such an environment.

Electromagnetic compatibility testing	IEC 60601 test level	Compliance level	Electromagnetic Environmental Guidelines.
Electrostatic discharge (ESD) according to IEC 61000 – 4 – 2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wooden or concrete or covered with ceramic tiles. If the floor is made of a synthetic material, the relative humidity must be at least 30 %.
Fast transient electric interference (burst) according to IEC 61000 - 4 – 4	± 2 kV for power lines	± 2 kV for power lines	The quality of the power supply voltage must correspond to a typical hospital or business power supply.
Voltage drops, temporary power outages, and variations in the power supply voltage according to IEC 61000 – 4 – 11	0 % reduction during 0,5 cycle 0 % reduction during 1 cycle 70 % reduction during 25 cycles 0 % during 250 cycles	0 % reduction for 2ms (drop) 0 % reduction for 4 ms (drop) 70 % reduction for 500 ms (drop) > 95 % reduction for 10 ms (drop)	The quality of the power supply voltage must correspond to a typical hospital or business power supply. If the user requires continued functionality even after failure of the power supply, we recommend operating the aScope™ Duodeno from an uninterruptible power supply.
Magnetic field at power supply frequency (50/60 Hz) according to IEC 61000-4-11	30 A/m (50 Hz)	30 A/m (50 Hz)	

Table 3 Electromagnetic Compatibility 2

Guidelines and manufacturer's declaration – resistance to electromagnetic interference. The aScope™ Duodeno is intended for operation in the electromagnetic environment specified below. The user of the aScope™ Duodeno should ensure that it is used in such an environment.

Electromagnetic compatibility testing	IEC 60601 test level	Compliance level	Electromagnetic Environmental Guidelines.
Line-conducted HF coupling according to IEC 61000 – 4 – 6	3 V;AM/1 kHz/ 80 % 150 KHz - 80 MHz	3 V	Portable and mobile radio equipment should not be used closer to the unit,
Electromagnetic fields according to IEC 61000 – 4 – 3	3 V/m; 80 MHz – 2.7 GHz 80 % 150 KHz - 80 MHz	3 V/m	including cables, than the recommended safety distance calculated according to the formula appropriate to the transmission frequency:
$d=3.5/3 \sqrt{P}$ up to 80 MHz $d=3.5/3 \sqrt{P}$ 80 MHz – 800 MHz $d=3.5/3 \sqrt{P}$ 800 MHz – 2.5 GHz where P is the nominal power of the transmitter in watts and d is the safety distance in meters.			

The field strength of stationary transmitters should be lower in all frequencies than the compliance level, according to an examination on site. Malfunctions are possible in the vicinity of equipment with the following symbols.

Table 4 Electromagnetic Transmission

Guidelines and manufacturer's declaration – resistance to electromagnetic interference. The aScope™ Duodeno is intended for use in the electromagnetic environment specified below. The user of the aScope™ Duodeno should ensure that it is used in such an environment.		
Transmission measurement	Compliance	Electromagnetic Environmental Guidelines.
HF transmission according to CISPR 11	Group 1	The aScope™ Duodeno uses HF power internally only. Its HF transmission is therefore very low and it is improbable that it will cause interference with electrical equipment in its vicinity.
HF transmission according to CISPR 11	Class A	The aScope™ Duodeno is suitable for use in facilities other than residential areas that are connected directly to the public power grid that also supplies buildings used for residential purpose, provided that the following warning is observed:
Upper harmonics according to IEC 61000 – 2 – 3	Class A	
Voltage fluctuations / flicker according to IEC 61000 – 3 – 3	Satisfied	
		<p>Warning:</p> <p>This device is intended only for use by trained medical personnel. This is a Class A device according to CISPR 11. In a residential area, this unit can cause radio interference, so it is necessary in this case to take suitable corrective measures, such as readjusting it, rearranging it, or shielding the unit or filtering its power connection.</p>


13. Contact

Manufacturer
 Ambu A/S
 Baltorpbakken 13
 2750 Ballerup, Denmark


Appendix 1. Connection Configurations Descriptions

The aBox™ Duodeno is equipped with multiple video outputs for viewing by clinicians and several options for documentation via additional video outputs and a remote signal (i.e. "trigger") to activate common recording devices. The steps in this appendix will guide the installer through various available connections and the required cabling for each configuration.


**Connecting the Mandatory Main Primary Imaging Monitor using the HD-SDI 1 Video Signal
Connecting the Optional Secondary Imaging Monitor using the HD-SDI 2 or DVI Video Signal**

Connecting of an Image Capture / Report Writing Workstation using the HD-SDI 2 or DVI Video Signal and TRIG 3

Connecting of a Digital Video Recorder using the HD-SDI 2 or the DVI Video Signal and TRIG 1 and/or TRIG 2

Appendix 2. WPF-Mediakit

WPF-Mediakit is an open source library.

It provides control for Visual Studio .NET to display video from a Windows imaging device.

<https://github.com/Sascha-L/WPF-MediaKit/wiki>

Version: 2.2.0

Release Date: 2017-01-19

Microsoft Public License (Ms-PL)

This license governs use of the accompanying software. If you use the software, you accept this license. If you do not accept the license, do not use the software.

1. Definitions

- The terms "reproduce," "reproduction," "derivative works," and "distribution" have the same meaning here as under U.S. copyright law.
- A "contribution" is the original software, or any additions or changes to the software.
- A "contributor" is any person that distributes its contribution under this license.
- "Licensed patents" are a contributor's patent claims that read directly on its contribution.

2. Grant of Rights

- **(A) Copyright Grant-** Subject to the terms of this license, including the license conditions and limitations in section 3, each contributor grants you a non-exclusive, worldwide, royalty-free copyright license to reproduce its contribution, prepare derivative works of its contribution, and distribute its contribution or any derivative works that you create.
- **(B) Patent Grant-** Subject to the terms of this license, including the license conditions and limitations in section 3, each contributor grants you a non-exclusive, worldwide, royalty-free license under its licensed patents to make, have made, use, sell, offer for sale, import, and/or otherwise dispose of its contribution in the software or derivative works of the contribution in the software.

3. Conditions and Limitations

- **(A) No Trademark License-** This license does not grant you rights to use any contributors' name, logo, or trademarks.
- **(B) If you bring a patent claim against any contributor over patents that you claim are infringed by the software, your patent license from such contributor to the software ends automatically.**
- **(C) If you distribute any portion of the software, you must retain all copyright, patent, trademark, and attribution notices that are present in the software.**
- **(D) If you distribute any portion of the software in source code form, you may do so only under this license by including a complete copy of this license with your distribution. If you distribute any portion of the software in compiled or object code form, you may only do so under a license that complies with this license.**
- **(E) The software is licensed "as-is." You bear the risk of using it. The contributors give no express warranties, guarantees or conditions. You may have additional consumer rights under your local laws which this license cannot change. To the extent permitted under your local laws, the contributors exclude the implied warranties of merchantability, fitness for a particular purpose and non-infringement.**

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