INSTRUCTIONS FOR USE

Single-use

- Ambu® Subdermal Corkscrew
- Ambu® Subdermal Twisted Pair Corkscrew
- Ambu® Subdermal Hooked Needle
- Ambu® Subdermal Curved Needle

Ambu

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Instructions for use

Recommended Application

Electrodes for neurophysiological examinations and intraoperative monitoring are intended for use with recording/monitoring equipment for the recording of biopotential signals including electroencephalogram (EEG), electromyograph (EMG) and nerve potential signals, and are intended for stimulation/recording with stimulation/recording equipment for electromyograph (EMG) and nerve potential signals.

Precautions

- Sterility is guaranteed up to the expiration date unless the package is opened or damaged. Do not use if individual pouch is damaged or has been opened.
- For single patient use only. Use on other patients can cause cross infection. Do not soak, rinse, or sterilize this device as these procedures may leave harmful residues or cause malfunction of the device. The design and material used are not compatible with conventional cleaning and sterilization procedures.

Warnings

- Do not connect directly to power source.
- Not to be used as ground electrodes.
- The electrode should be used by health care professionals trained in performing neurophysiological examinations and intraoperative monitoring - US: Rx only.
- The electrode is for use in hospital/clinic settings.
- During procedures where electro surgery is used, attention should be drawn to the risk of unintended interference with the electrode, which could cause unintended heating of the needle and damage to the tissue. To avoid this the electrodes should be removed or placed as far as possible from the electro-surgical area.
- If used as stimulation electrodes, thermal burns could result with prolonged stimulation outputs above 0.25W/cm².
- If used as recording electrodes, the possible loss of electrocautery return path may result in high current density passing through these needles with possible tissue heating and damage resulting.

Directions for use

Electrode is compatible with IEC 60601-1 compliant Electroencephalography (EEG), electromyograph (EMG), or intraoperative monitoring (IOM) monitoring equipment with DIN 42-802 input socket.

When selecting a disposable needle electrode, it is important to use one that is of suitable shape and diameter and wire length for the investigation being undertaken.

1. Peel the sterile pouch open and remove the protection cap/tube if applicable.
2. Connect the needle electrode connector to the monitoring device. Verify that the electrode is connected correctly to the monitoring device.
3. The correct placement of needle electrodes is typically determined by the hospital/medical/and/or local protocols/procedures.
4. Apply the electrode by placing it against the scalp and turning clockwise 1/2 to 1 complete turn (only for corkscrew).
5. Remove by turning the electrode counterclockwise (only for corkscrew).
6. After use, dispose the needle electrode in a designated container.

US: Rx only

For symbol explanation see www.ambu.com or www.ambusa.com

Ambu® is a registered trademark of Ambu A/S, Denmark
Ambu A/S is certified according to ISO 13485.

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