Directions for Use

Ambu® CPR-Software

Also
designed
for WLAN

Contend

1	Introduction	3
2	Getting Started	3
2.1	System Requirements (PC)	3
2.2	System Compatibility (Ambu equipment)	4
2.3	Installation of the CPR-Software	4
2.3.1	Windows Vista, 7, 8	4
2.3.2	Windows XP	4
2.4	Installation of the USB Driver for the Ambu manikins	5
2.4.1	Windows Vista, 7, 8	5
2.4.2	Windows XP	6
2.5	Power supply of Ambu manikins	7
3	Connections	7
3.1	Single Manikin – USB	7
3.2	Single Manikin – WLAN	7
3.3	Single Manikin – model C manikin	8
3.4	Multi Manikin – USB	9
3.5	Multi Manikin – WLAN	9
3.6	Multi Manikin – model C manikins	10
4	CPR Software functionalities	10
4.1	Configuration window	10
4.2	New session window	11
4.3	Add Manikin	12
4.3.1	USB manikins	12
4.3.2	WLAN manikins	13
4.4	Main Window	13
4.5	The graph window	15
4.6	Debriefing	16
4.7	Course planer	18
4.8	Surveillance Window	19
4.9	CPR Software Database	20
4.10	ECG-Box	21
4.11	Getting Help	23
5	Configuration of Ambu manikins	23
5.1	Configuration of the WLAN	24
5.2	Default settings	24
5.3	Reset to default	24
5.4	Firmware update	25
6	Troubleshooting	25

1 Introduction

The CPR software is the ideal tool to monitor and register your training with Ambu manikins. It provides you with all the essential training parameters so that you can control and improve the training sessions. It is the possibility of saving and printing all information that helps you grade a trainee's performance and later to document your decision.

It combines all features of the well-known and accepted former versions of the software and in addition makes it possible to connect to the Ambu manikins directly through USB or using a wireless local area network (WLAN) connection.

The software is designed with the idea of a simple and structured display where a trainer can survey and record the performance of one or multiple teams on one PC. The single sessions are stored so that they can be evaluated afterwards.

The Ambu CPR-Software can be used with all Ambu manikins that are equipped with a computer connection.

2 Getting Started

2.1 System Requirements (PC)

Minimum computer requirements

Computer: Pentium Processor

RAM: 8 MB

Operating System: Windows XP, Vista, 7, 8

Free space on Hard Drive: 10 MB

Screen resolution: 640x480, 256 color.

1 free USB-port

1 free COM-port (for model C manikins)

Miscellaneous Mouse or pointing device

Keyboard

Recommended computer requirements

Computer: Pentium Processor

RAM: 16 MB

Operating System: Windows XP or Windows 7

Free space on Hard Drive: 20 MB

Screen resolution: 1024x768, 32 bit color.

1 free USB-port

WLAN Adapter

Mouse or pointing device

Keyboard

2.2 System Compatibility (Ambu equipment)

The Ambu CPR Software is designed to be connected directly to one of the following Ambu manikins:

- Ambu® AmbuMan W
- Ambu® Defib Trainer System W
- Ambu® Airway Man W
- Ambu® Cardiac Care Trainer System W

The Ambu CPR-Software is also backward compatible to the manikins produced before 2009 known as the models C.

To connect the model C manikins to the Ambu CPR-Software the Ambu CPR-Software Kit is required. The Ambu CPR-Software Kit is available in the following configurations:

- Ambu CPR Software Kit (for RS232 connection)
- Ambu CPR Software Kit with USB Interface

2.3 Installation of the CPR-Software

Before you install the Ambu-CPR Software make sure you have all administrative rights for the installation. The Ambu CPR-Software will need access to the ODCB-Database from Microsoft® Windows.

2.3.1 Windows Vista, 7, 8

The Ambu CPR Software is delivered on a USB drive and installs like most other Windows applications.

To install the Ambu CPR Software from the USB drive, do the following:

1. If Windows is not running, please start it now.
2. Insert the Ambu CPR Software USB drive.

Windows should detect the USB drive and automatically run the installation software. If the automatic detection feature is turned off, please start the installation software by clicking Start and then on "Computer". Select the USB Drive and double-click on the file "setup_cpr.exe".

3. Click on "Allow" if Windows shows the screen "An unidentified program wants to access to your computer".
4. Follow the instructions on the screen.

During the installation you have the possibility of choosing the language version you would like to install. It is possible to mark and hereby to install more than one language.

2.3.2 Windows XP

The Ambu CPR Software is delivered on a USB drive and installs like most other Windows applications.

To install the Ambu CPR Software from the USB drive, do the following:

1. If Windows is not running, please start it now.

2. Insert the Ambu CPR Software USB drive.

Windows should detect the USB drive and automatically run the installation software. If the automatic detection feature is turned off, please start the installation software by clicking Start and then on "Computer". Select the USB Drive and double-click on the file "setup_cpr.exe".

3. Follow the instructions on the screen.

During the installation you have the possibility of choosing the language version you would like to install. It is possible to mark and hereby to install more than one language.

2.4 Installation of the USB Driver for the Ambu manikins

Before you install the USB-Drivers for the Ambu manikins make sure you have all administrative rights for the installation.

NOTE: If you have more than one USB port and you are using the manikins on different USB ports, the installation of the USB driver for Ambu manikins has to be repeated on each USB port.

2.4.1 Windows Vista, 7, 8

To install the USB drivers for the Ambu manikins do the following:

1. Insert the Ambu CPR Software USB drive.
2. Connect a manikin to the USB port and turn the manikin on by pressing the power button.
3. Windows will detect the manikin and will automatically ask you to install the driver.
4. Click on "Locate and install driver software (recommended)".
5. Click on "continue" if Windows shows the screen "Windows needs your permission to continue".
6. Click on "install this driver software anyway" if Windows shows the screen "Windows can't verify the publisher of this software"
7. To finish the installation click on "close"

If you want to install the Ambu USB Interface for Ambu manikins model C do the following:

Note: Windows will install a driver for the Ambu Controller and a driver for Ambu AMC Box. Both drivers must be installed for the USB Interface to work.

1. Insert the Ambu CPR Software USB drive.
2. Connect the USB Interface to the USB cable and to the computer.
3. Windows will display the screen "Found new Hardware - Ambu USB Interface" and will ask for a driver. (If your computer does not find the USB Interface you have disabled the automatic installation for USB hardware. Please enable this Windows function).
4. Click on "Locate and install driver software (recommended)".

5. Click on "continue" if Windows shows the screen "Windows needs you permission to continue".
6. To finish the installation click on "close".
7. Windows will now display the screen "Found new Hardware – USB Serial Port" and will ask for a driver.
8. Click on "Locate and install driver software (recommended)".
9. Click on "continue" if Windows shows the screen "Windows needs you permission to continue".
10. To finish the installation click on "close".

2.4.2 Windows XP

To install the USB drivers for the Ambu manikins do the following:

1. Insert the Ambu CPR Software USB drive.
2. Connect a manikin to the USB port and turn the manikin on by pressing the power button.
3. Windows will detect the manikin as and will automatic ask you to install the driver. (If your computer does not find the USB Interface your have disabled the automatic installation for USB hardware. Please enable this Windows function.)
4. Click on "No, not this time" if Windows asks to connect to Windows Update to search for software and click on "next".
5. Click on "Install the software automatically (recommended)" and click on "next".
6. Click on "continue anyway" if Windows shows a screen saying the software has not passed Windows Logo testing.
7. To finish the installation click on "Finish".

If you want to install the Ambu USB Interface for Ambu manikins model C do the following:

Note: Windows will install a driver for the Ambu Controller and a driver for Ambu AMC Box. Both drivers must be installed for the USB Interface to work.

1. Insert the Ambu CPR Software USB drive.
2. Connect the USB Interface to the USB cable and to the computer.
3. Windows detect the Ambu Controller and will display the screen "Found new Hardware" and will ask for a driver. (If your computer does not find the USB Interface your have disabled the automatic installation for USB hardware. Please enable this Windows function).
4. Click on "No, not this time" if Windows asks to connect to Windows Update to search for software and click on "next".
5. Click on "Install the software automatically (recommended)" and click on "next".
6. Click on "continue anyway" if Windows shows a screen saying the software has not passed Windows Logo testing.
7. To finish the installation click on "Finish".
8. Windows will now detect the Ambu AMC Box and will display the screen "Found new Hardware" and will ask for a driver. (If your computer

- does not find the USB Interface you have disabled the automatic installation for USB hardware. Please enable this Windows function).
9. Click on "No, not this time" if Windows asks to connect to Windows Update to search for software and click on "next".
 10. Click on "Install the software automatically (recommended)" and click on "next".
 11. Click on "continue anyway" if Windows shows a screen saying the software has not passed Windows Logo testing.
 12. To finish the installation click on "Finish".

2.5 Power supply of Ambu manikins

Ambu manikins (model W) need a power source to be used with the CPR Software. The power supply can be established by using the USB cable, batteries or the optional AC/DC Adapter.

For recommendations and instructions please look at the direction for use of the manikin.

3 Connections

To connect an Ambu manikin to the computer you must have the Ambu CPR Software and the USB Drivers installed. The connection to the manikin can be made by USB or WLAN.

If you are working with an Ambu manikin model C you have to use the ADC Cable or the Ambu USB interface that comes with the CPR-Software Kit.

3.1 Single Manikin – USB

Release the monitoring instrument on the manikin and attach the USB cable to the USB socket on the left side of the monitoring instrument. Connect the other end of the USB cable to the Computer. The yellow LED above the USB socket on the monitoring instrument will indicate when the connection is established.

Turn on the manikin by pressing the power button on the right side of the monitoring instrument. The green LED above the button will start flashing while the system is booting up. Once the system is ready the LED will flash steady.

3.2 Single Manikin – WLAN

The Ambu manikins (model W) are using a standard WLAN module and are preconfigured to be running an ad-hoc network. The settings of the WLAN can be changed as described in 5.1. We recommend that the settings be only changed if you have good knowledge about computer networking.

Release the monitoring instrument on the manikin and turn on the manikin by pressing the power button on the right side of the monitoring instrument. The green LED above the button will start flashing while the system is booting up. Once the system is ready the LED will flash steady.

Turn on the WLAN adapter on your computer and open the Windows Wireless network connection. Search for wireless networks. The default network name for Ambu manikin is "Ambuwireless". Click on connect to connect to the manikin. Windows will now establish the wireless connection.

NOTE: If a wireless connection to a manikin is chosen, make sure the network setting (like network name and encryption / password) for the manikins are correct. Wrong settings can prevent the computer to be able to connect to the manikin.

Use the Windows wireless connection manager to connect to the same network as the manikin is connected to. If the manikin and the PC are connected to different networks the CPR software will not be able to detect the manikin.

Windows 8.1

If the standard ad hoc connection should be used, click on "create Ambu ad hoc network" on your App screen. Follow the instruction.

To connect to the new created network, click on "Connect to Ambu ad hoc network". Windows will now establish the wireless connection.

3.3 Single Manikin – model C manikin

The Ambu manikin models C can be connected using the ADC cable or the USB interface.

If you are using the ADC cable connect the 6-pin DIN plug to the connector on the plug on the left side of the monitoring instrument on the manikin. Connect the 9-pin plug to a serial RS232 port of the computer. If necessary you can place an RS232 extension of max. 10 m between the computer and the ADC cable.

If you are using the Ambu USB Interface connect the 6-pin DIN plug to the connector on the plug on the left side of the monitoring instrument on the manikin. Connect the 25-pin plug to the USB interface. Then

connect the USB cable to the USB Interface and the other end of the USB cable to a USB Port on the Computer.

3.4 Multi Manikin – USB

Multiple manikins can be used together with the CPR Software. That maximum number of manikins used at the same time should not exceed 9 manikins.

For every manikin a free USB port is needed. If your computer does not have enough free USB ports an external USB HUB with external power supply can be used. The USB Hub must be capable to deliver constantly 500 mA to each USB Port.

Note: The USB Hub must be capable to deliver constantly 500 mA to each USB Port. The power supply must therefore be high enough to deliver the current for the USB port and the Hub itself.

To connect the manikins do the following:

Release the monitoring instrument on the manikin and attach the USB cable to the USB socket on the left side of the monitoring instrument. Connect the other end of the USB cable to the Computer or an external USB HUB. The yellow LED above the USB socket on the monitoring instrument will indicate when the connection is established.

Turn on the manikin by pressing the power button on the right side of the monitoring instrument. The green LED above the button will start flashing while the system is booting up. Once the system is ready the LED will flash steady.

Repeat this procedure for all manikins that shall be used.

3.5 Multi Manikin – WLAN

Multiple manikins can be used together with the CPR Software. That maximum number of manikins used at the same time should not exceed 9 manikins.

The Ambu manikins (model W) are using a standard WLAN module and are preconfigured to be running an ad-hoc network. The settings of the WLAN can be changed as described in 5.1. We recommend that the settings be only changed if you have good knowledge about computer networking.

Release the monitoring instrument on the manikin and turn on the manikin by pressing the power button on the right side of the monitoring instrument. The green LED above the button will start flashing while the system is booting up. Once the system is ready the LED will flash steady.

Repeat this procedure for all manikins that shall be used.

Turn on the WLAN adapter on your computer and open the Windows Wireless network connection. Search for wireless networks. The default network name for Ambu manikin is "Ambuwireless". Click on connect to connect to the manikins. Windows will now establish the wireless connection

NOTE: If a wireless connection to a manikin is chosen, make sure the network setting (like network name and encryption / password) for the manikins are correct. Wrong settings can prevent the computer to be able to connect to the manikin.

Use the Windows wireless connection manager to connect to the same network as the manikin is connected to. If the manikin and the PC are connected to different networks the CPR software will not be able to detect the manikin.

Windows 8.1

If the standard ad hoc connection should be used, click on "create Ambu ad hoc network" on your App screen. Follow the instruction.

To connect to the new created network, click on "Connect to Ambu ad hoc network". Windows will now establish the wireless connection.

3.6 Multi Manikin – model C manikins

Multiple manikins can be used together with the CPR Software. That maximum number of manikins used at the same time should not exceed 9 manikins.

For every manikin a free USB port is needed if the Ambu USB interface is used. If your computer does not have enough free USB ports an external USB HUB with external power supply can be used.

For every manikin a free serial RS232 port is needed if the ADC Cable is used. If your computer does not have enough free serial RS232 ports a serial port extension can be used.

If you are using the ADC cable connect the 6-pin DIN plug to the connector on the plug on the left side of the monitoring instrument on the manikin. Connect the 9-pin plug to a serial RS232 port of the

computer. If necessary you can place an RS232 extension of max. 10 m between the computer and the ADC cable.

If you are using the Ambu USB Interface connect the 6-pin DIN plug to the connector on the plug on the left side of the monitoring instrument on the manikin. Connect the 25-pin plug to the USB interface. Then connect the USB cable to the USB Interface and the other end of the USB cable to a USB Port on the Computer.

Repeat this procedure for all manikins that shall be used.

4 CPR Software functionalities

You can start the Ambu CPR Software from the Start menu

- Click the "Start" button. The Start menu opens.
- Choose "Programs". The Programs folder opens.
- Choose "Ambu CPR Software". The Ambu CPR Software folder opens.
- Click "Ambu CPR Software" (desired language)

4.1 Configuration window

To change the configuration of the Ambu CPR Software click on "Tools" in the menu bar and then click on "Configuration".

The configuration windows is separated into 4 Tabs:

- General
- Debriefing 1
- Debriefing 2
- Print
- CPR Controller

To change settings click on the tab on the top of the window and do the necessary changes in the single fields.

In the "General" setting you can configure if you want to play a sound when the session is finished. You can set the frequency of the metronome that can be activated during the session and if you generally want to show or hide the ECG-signal in the session.

Note: The ECG signal can only be displayed if a Defib Trainer W or a Cardiac Care System W is connected.

The tab "Debriefing 1" shows you the settings of the ventilation parameters and how they are displayed in the debriefing view.

The tab "Debriefing 2" shows you the settings of all compression parameters and how they are displayed in the debriefing view.

In the tab "Print" you can change the setting for the printout of the sessions and you can insert a text for the page heading.

If you want to print the graph you can set the maximum time of the graph that will be printed per page.

If you are using the Ambu CPR Printer you adjust and configure the CPR controller by clicking on the tab "CPR Controller". In order to transfer the data to the CPR Controller it must be connected to the computer. For further detail please see the directions for use for the Ambu CPR Printer.

4.2 New session window

When you start the Ambu CPR Software, a routine starts to detect all manikins that are connected with wires to the computer. If you have connected model C manikins the CPR- Software will automatically calibrate the manikins. Model W manikins do not need to be calibrated as they have an build-in automatic calibration function insider the manikins.

A Scenario consists of an Algorithm, the Duration and which Training mode you want to use in a session.

The pull-down menu gives you the possibility to choose a predefined scenario. It is possible to modify the chosen scenario with the Algorithm, Duration, and Training mode pull-down menus. By clicking the "Details" buttons, you can create or delete an existing scenario or algorithm.

The software automatically detects all wired connected manikins. You can choose the manikins by using the drop down menu on the screen.

The number of participants varies depending on the chosen training mode.

The software provides you with the possibility to personalize and organize the training session. Out of a database that you access via the Details buttons it is possible to choose the participants for the Ventilation and/or Compression training.

Press "Continue" to enter the Main Window

4.3 Add Manikin

Click on “Add Manikin” if you want to add a manikin, you are using wireless connected manikin or you want to re-calibrate the connected model C manikins.

4.3.1 USB manikins

If you want to add a manikin on that is connected by USB or you have to calibrate a model C manikin click on “USB” to search for the manikins.

Model W manikins will automatically be detected and do not need to be calibrated.

Model C manikins will automatically be detected but need to be calibrated before use. Checkmark the manikins that need to be calibrated and click on “Calibrate selected manikins”.

Click on “Close” to return to the “New session” window.

4.3.2 WLAN Manikins


If you want to add a manikin that is connected by WLAN click on “WLAN / LAN” to search for the manikins.

The manikins will automatically be detected and do not need to be calibrated.

The default name of the WLAN manikins is "Ambu_Default". You can change the name as described in 5.1

Click on "Close" to return to the "New session" window.

4.4 Main Window

4.5 The graph window

The graph window is used to visualize the different training parameters (A-H). Furthermore it provides you with on-screen evaluating features (I,J) and additional information (K-P)

Ventilation features/information

- A) Ventilation volume
- B) Stomach inflation indicator

Compression

- C) Actual compression Rate (1/min)
- D) Compr./Relax. Ratio
- E) Compression depth
- F) Wrong hand position
- G) Leaning
- H) Cardiac output

Event management

- I) Event toolbar
- J) Event trace (not shown)

Evaluating

- K) Modify displayed time interval (Zoom in/out)
- L) Compression-ventilation coordinates
- M) Actual recommendations
- N) Actual manikin
- O) Name of the participants
- P) Start/pause a session
- Q) Actual time out of total time
- R) Link to the result window (Debriefing)
- S) ECG signal
- T) show/hide ECG Box
- U) Battery status

It is possible to measure time intervals by clicking on the left mouse button and dragging the area that should be evaluated


By clicking the right mouse button you can decrease the time display. For further information please refer to the information in the online-help (see chapter 4.11).

Pause the session and press "Debriefing" to enter the Debriefing view.

Note: The ECG-Signal and the ECG-Box are only available if you connect an Ambu Defib Trainer W or an Ambu Cardiac Care Trainer System W.

4.6 Debriefing

The debriefing view makes it possible for the trainer to form a general view of the students ability to carry out CPR.


The indicator (A) must be within the green section. This indicates that the various parameters have been carried out correctly according to the recommendations in force.

(B) The values shown in this column display the average or the percentage achieved during the training session.

If you want to add personal recommendation you can use the text field (C).

For further information please refer to the information in the online-help (see chapter 4.11).

Press the **Review** button to proceed to the final step of a session.


When the session has been prepared, visualized, and the debriefing is over, the final step of the session process is the review view.

In the review view, it is possible to change the names of the participants in the session, to mark whether the metronome was used or not, to decide whether the session is approved or not, and to mark who approved the session. Finally, it is possible to write notes regarding the session.

When all of this has been completed, you should press "Close" (if you paused a session the button is named "Abort") to mark that the session is over. The software will now ask you if you want to save the session in the database.

4.7 Course planer


When you connect many Ambu manikins to the computer, we recommend that you use the Course Planner to plan the training sessions.

When you plan a course, first you should enter your name as the instructor of the course. Then, select which scenario should be the default one for this group session. Now, for each training session, enter the names of the participants in the next empty row, and select which manikin, the session is performed with. The software will automatically fill in preliminary start and end-times of the session, and prepare a new, empty row for the next training session.

If one of more of the training sessions should be done with other parameters than the default scenario, press the Details button in the corresponding row, and change the parameters as wanted.

When all of the course has been planned, press the "Run Course" button to proceed. If you want to, you can press "Print" to get an overview of the session printed on your printer.

In the Run course window, the software is ready to handle the training session.

Advertise the participants of each session, and when they are ready, press the "Play" button in the corresponding row to start the training session. A graph view will appear behind the Run course window, and register the training session in the background in the graph view.

When a training session completes, a check mark will appear in the corresponding row of the Run course window, and you can prepare for the next training session.


During the execution of a course, you can make last minute changes in the order of sessions or setup of each session with the Move row up, Move row down, and Details buttons.

By using the toolbar icons, you can move back and forth between the graph view or the Run course window and monitor each individual session in detail, or the overall progress of the course. You can also use keyboard short-cuts to speed up this navigation. See the online help for more information.

In addition to these facilities, the software provides a surveillance window which provides an overview of the graphs of many training sessions at once.

4.8 Surveillance Window

This window allows you to survey the different sessions, you are running.


Active sessions are shown on a white background, while completed sessions have a grey background.

By double clicking database dialogs button, or single clicking the right mouse button, on one of the views, you open the corresponding graph window.

4.9 CPR Software Database

The database dialogs in the program are designed with a standard interface. They share the same appearance and have the same buttons for navigation and editing.

For example, the **Database | Algorithms** dialog:

At the top of the dialog, you will find the Quick Find drop down list. Use this list to quickly locate any given algorithm by clicking on it, scroll in the list until you find the one you need to inspect, or edit. This principle applies to all database dialogs. The following directions will apply to the all of the persons, scenarios, and courses dialogs.

In order to understand how the database system works, imagine that all of the algorithms are lined up next to each other in an alphabetically ordered list, where one of these is the current record shown on screen.

When you use the four navigation buttons next to the Quick Find list, you can move to the previous record, the next record, the first record or the final record. This record now becomes the current record.

When you change fields in the algorithm, you only edit the current record. If you change the name of the algorithm, you will change its position in the list (this can not be seen on screen, because the algorithm is still current):

We don't recommend you to change the pre-defined algorithms, so press "Discard changes" to drop the changes, and the record will return to the original position in the original form.

If you press "New algorithm", a new algorithm will be created at the extreme right of the list. This becomes the active algorithm.

4.10 ECG Box

The ECG Box is only available if you connect an Ambu Defib Trainer W or an Ambu Cardiac Care Trainer System W.

There are 26 different ECGs.

The upper row simulates the 6 basic ECGs which are especially important in training emergency situations:

The four buttons surrounded by a white box are the control buttons. These buttons are used by the instructor to manage the sequences of the different cardiac rhythms. The ECG rhythm currently generated in the manikin is indicated by a constant green light diode. When clicking on another ECG button the corresponding diode will start to flash and this rhythm is now on standby.

CHANGE RHYTHM:

By clicking this button the instructor can change the current ECG rhythm to the ECG rhythm on standby. The former ECG rhythm switches off completely.

SWOP RHYTHMS:

The instructor can swap between the current ECG rhythm and the standby ECG rhythm. For example after activating VENTR. TACH the instructor can select NORMAL SINUS thus it will be on standby mode (the diode flashes). When clicking SWOP, the NORMAL SINUS will take over and the VENTR. TACH will be on standby mode. In this way the instructor can swap between two different cardiac rhythms. Accordingly, strokes of certain cardiac rhythms can be trained.

NO DEFIB. EFFECT:


When activating this button, defibrillation on the Ambu Cardiac Care Man does not have any effect, even if the instructor has selected another rhythm as standby. For example: VF COARSE or FINE is activated. NO DEFIB EFFECT is active and NORMAL SINUS is on standby. When defibrillation is carried out, VF continues, because NO DEFIB EFFECT is activated. When clicking NO DEFIB EFFECT again it is deactivated and the diode is turned off. The next time defibrillation is carried out the ECG switches to NORMAL SINUS.

PULSE ON/OFF:

When the button "Pulse on/off" is activated, the enabled pulse signal will generate a palpable pulse on the Ambu Cardiac Care Man except for the rhythms VF COARSE, VF FINE, ASYST. and IDIOV. RHYTHM.

The diode will flash in a pulse rhythm if the selected cardiac rhythm has a pulse. If the diode does not flash but is on constantly there is no pulse. If "Pulse on/off" is deactivated a palpable pulse will not be generated on the manikin and EMD or PEA can be simulated (EMD: Electro Mechanical Dissociation, PEA: Pulseless Electrical Activity).

The ECG Box can be shown/hidden by clicking on the ECG Box Button. It is possible to move the ECG Box around the screen.

In the surveillance view the ECG Box is linked to a specific manikin. You can choose the manikin you want to use the ECG Box for by clicking on the drop down menu "Show ECG Box for".

4.11 Getting Help

<F1> Button

You can get help regarding a specific feature by placing the cursor over the feature and pressing <F1>.

Toolips and Status Bar

To see the name and function of a button on the toolbar, or a feature in a palette, place the cursor over it. A short description appears after a short time.

Help Buttons in Dialog Buttons

If you have questions about a dialog box, click the Help button in the dialog box.

Online Help

Choose "Help Contents" from the Help menu to access the contents, index, and keyword search of the online Help.

Ambu Homepage

If you have Internet access, please visit our Ambu Web page on www.ambu.com. Here you will find answers to frequently asked questions, Software updates, and other interesting information.

5 Configuration of the manikins

For a personal use of the manikins the setting for the WLAN connection can be adjusted.

To change the settings connect the manikin through USB and start the configuration software "WiFi.exe" from the CD or the program folder.

5.1 Configuration of the WLAN

The following settings can be changed to adjust the manikin to the personal needs.

We recommend that only persons with good knowledge in computer networking do the configuration of the manikins.

5.2 Default settings

The default settings can be seen on the screen below.

Note: The default setting is an open and not secured wireless network.

5.3 Reset to default

If the manikin needs to be reseted to default values click on the button "Reset" in the configuration software or press the power button on the manikin for 10 seconds.

5.4 Firmware update

If a firmware update is necessary connect the manikin to a computer through USB. Start the firmware update file by double clicking on it. The update will be done automatically and you will get information when the update is finished.

The newest firmware version can be found at www.ambu.com.

6 Troubleshooting

Problem	Reason	Solution
The Ambu® CPR software cannot find the manikin	- USB driver is not installed	Please install the USB Driver
The "Ambu Controller" is not shown in the Device Manager under "Universal Serial Bus Controller"	- USB driver is not installed - USB port was changes	Please install the USB Driver Please connect the manikin to the USB Port where the Interface was originally installed
The Ambu® CPR software cannot find a manikin although I have installed the drivers	- Another USB driver has overwrite the driver for the manikin - USB port was changes	Please reinstall the USB Driver Please connect manikin to the USB Port where the Interface was originally installed
My Computer says that the USB hardware cannot be used	- Another USB driver has overwrite the driver for the manikin	Please reinstall the USB Driver

	<ul style="list-style-type: none"> - USB port was changes 	Please connect manikin to the USB Port where the Interface was originally installed Please connect the USB Interface to the USB Port where the Interface was originally installed Please reinstall the USB Driver
In the Device Manager under "Universal Serial Bus Controller" is a yellow mark in front of the "Ambu Controller"	<ul style="list-style-type: none"> - USB port was changes 	
The "Ambu Controller" with a yellow mark in front is shown under "Other Devices" in the Device Manager under "Universal Serial Bus Controller"	<ul style="list-style-type: none"> - Another USB driver has overwrite the driver for the USB Interface -Windows has deleted the drivers for the USB Interface 	Please reinstall the USB Driver
In the Device Manager under "Ports" I don't have an entry "Ambu AMC Box"	<ul style="list-style-type: none"> - Another USB driver has overwrite the driver for the USB Interface 	Please reinstall the USB Driver
The text on the ECG Box does not fit the buttons No wireless connection can be established	<ul style="list-style-type: none"> - The font is not installed correctly - Computer and manikin is in different networks - The settings of the manikins and the computer are not the same 	<p>Restart the computer after the installation.</p> <p>Open the wireless connection and search for networks. Connect the computer to the network of the manikin Set the manikin and the computer to the same settings.</p>

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Um die Ziffern anzeigen zu lassen, auf welchen die Auswertung basiert, drücken Sie auf die Tasten „Details“ ©.
Ambu A/S