

aCart™ Compact

DATASHEET

aCart™ Compact is a modern, mobile workstation that adds value to the Ambu single-use product portfolio. The aCart Compact workstation provides convenient storage and portability without taking up a lot of space.

A modern mobile workstation

The aCart Compact workstation enables you to move fast when you need to respond to an emergency or change clinical environments. It is also ideal for the outpatient clinic or for doctor's offices with limited space.

Robust and portable

The design and choice of materials make the aCart Compact workstation light enough to easily transport, yet robust and stable enough to protect your valuable endoscopic equipment. The hygiene of the aCart Compact workstation is certified by the DGKH (German Society for Hospital Hygiene).



KEY BENEFITS

- Height adjustable IV pole
- Ergonomic design handle/ cable winder
- Monitor/displaying and processing unit holder
- Preparation tray/shelve
- Storage hook for spare endoscopes
- Support column with integrated cable duct
- Storage basket for accessories, cables or airway tools
- Base with four lockable double swivel casters

Ambu

TECHNICAL SPECIFICATIONS

Dimensions [mm (in)]

Width	530 (20.86)
Height without IV pole	1320 (51.97)
Height with IV pole min. / max.	1511 (59.49) / 2351 (92.55)
Depth	610 (24.02)
Package Dimensions W/D/H	1180/590/225 (46.46/23.23/8.86)

Weight [kg (lbs)]

Total	23 (50.71)
-------	------------

Load capacity [kg (lbs)]

Maximum load	35 (77.1)
Shelve	10 (22)
Monitor holder	14 (30.8)
Basket	3 (6.6)
IV Pole	5 (11)
Optional Vario Mount	8 (17.6)

Device classification

Europe	Class I, CE certified according to DIN EN 60601-1
USA	Class I, 510(k) exempt from premarket notification

Operating environment

Operating temperature [°C (°F)]	10 - 40 (50 - 104)
Operating relative humidity [%]	30 - 75
Operating atmospheric pressure [kPa]	70 - 106

Storage environment


Storage temperature [°C (°F)]	-25 - 70 (-13 - 158)
Storage relative humidity [%]	10 - 95
Storage atmospheric pressure [kPa]	50 - 120

ORDERING SPECIFICATIONS


Description	Reference number
aCart™ Compact	KU.9713.800

Spare parts & Accessories

aCart™ Compact Vario Mount	HA.3000.991
aCart™ Compact Basket	KN.1244.991
aCart™ Compact Universal Endoscope Holder	KN.1245.800
aCart™ Compact Drawer	KN.1171.903
aCart™ Compact Shelf	KU.9716.991



aCart™ Compact (Right)



aCart™ Compact (Left)

Ambu

Distributor:

Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
ambu.com

Manufacturer:

iTD GmbH
Jahnstrasse 1
84347 Pfarrkirchen
Germany